

Medical device product instructions

Medical ultrasound coupling agent

[Product Name] Medical Ultrasound Coupling Agent

[Models and Specifications] 5g, 10g, 15g, 20g, 25g, 30g, 50g, 60g, 80g, 90g, 100g, 110g, 120g, 150g, 160g, 180g, 200g, 250g, 260g, 180g, 300g, 400g, 500g, 800g, 1000g, 2000g. Special specifications can be customized upon request.

[Product Description] Composed of carbomer, glycerin, triethanolamine, and water. During ultrasound diagnostic or therapeutic procedures, this product is used as an intermediary medium between the skin and the radiating surface of the probe (or treatment head) to transmit sound waves. It is not sterile and is not sterile.

[Intended Use] Improves ultrasonic coupling between the probe and the patient. For use on intact skin.

[Labeling and Marking Instructions]

			
Avoid sun exposure	Note, refer to the accompanying documents	Keep dry	Stacking layer limit

[Instructions] For external use. Open the package and apply to intact skin.

[Precautions, Warnings, and Tips]

1. Do not use if the packaging is damaged.
2. Do not use if leakage, mold, turbidity, or other signs of deterioration are observed.
3. Keep out of reach of children. Do not take orally.
4. Do not use if you are allergic to the product.

[Storage] Store in a dry place at room temperature.

[Expiry Date] Three years.

[Product Registration Number/Product Technical Requirements Number]

Xiangchang Medical Device Registration Number 20240011

[Production Registration Number] Xiangchang Medical Device Production Registration Number 20210037

[Manufacturing Batch Number, Production Date, Expiration Date] See packaging

[Registered Company/Manufacturer] Changde Yinuo Biopharmaceutical Co., Ltd.

[Registered Address/Manufacturing Address] Fengkang Avenue, Xicheng Community, Dajinggang Town, Anxiang County, Changde City, Hunan Province (within Jinggang Park)

[After-Sales Service] Changde Yinuo Biopharmaceutical Co., Ltd.

[Date of Instruction Sheet Preparation] July 11, 2024